



Test Gri3D™

Grow organoids reproducibly by the hundreds


From single cells



From aggregates of cells



Continuous growth analysis from time of cell seeding



Microwell vs. seeding density guide

		Microwell diameter						
		10	50	100	300	500	1000	1500
Starting cell number per aggregate	single cell	✓	✓	✓	✓	✓	✗	✗
	2-20	✗	✓	✓	✓	✓	✗	✗
	20-50	✗	✗	✓	✓	✓	✗	✗
	50-100	✗	✗	✗	✓	✓	✓	✗
	100-500	✗	✗	✗	✗	✓	✓	✗
	500-2000	✗	✗	✗	✗	✗	✓	✓
	3000+	✗	✗	✗	✗	✗	✗	✓

Number of microwells per 24 well plate well

		Microwell diameter						
		10	50	100	300	500	1000	1500
Gri3D™ size	ø 6mm	80'000	6'500	2'200	300	100	19	7
	ø 8mm	145'000	11'500	4'000	550	200	37	19
	ø 10mm	225'000	18'500	6'200	880	300	80	37